

DECLARATION OF CONFORMITY

Manufacturers : HUBIT Co., Ltd.
(Ojeon-dong, Byeoksantechnopia), A-1204, 13 Ojeongongeoop-gil,
Uiwang-si, Gyeonggi-do, Korea

EC Representative : SC ORTOFORUM SRL
224 ANDREI SAGUNA STREET, 310077 ARAD,
ROMANIA

Product Name : Orthodontic Brackets

Brand Name : PERFECT CLEAR, PERFECT CLEAR II,
PERFECT CLEAR II (TRI-EDGE) and WOW Bracket
(Refer to the attachment for detailed type numbers.)

Classification : Class IIa (according to MDD Annex IX, Rule 5)

GMDN Code : 31759

The product identified above complies with the essential requirements of the above EC Directives by meeting the following standards:

Standards applied : EN ISO 15223-1:2008, EN 1041:2008, EN 1641:2009, ISO 7405:2008,
EN ISO 10993-1:2009/AC:2010, EN ISO 10993-5:2009, EN ISO 10993-10:2010, EN ISO 10993-11:2009,
EN ISO 13485:2012/AC:2012, EN ISO 14155:2011, EN ISO 14971:2012, MEDDEV 2.7.1 Rev.3, MEDDEV
2.12.1 Rev.7, MEDDEV 2.12.2 Rev.2

*This Declaration of Conformity is based on the EC Directives 93/42/EEC amended by 98/79/EC, 2000/70/EC, 2001/104/EC, 1882/2003, Annex II (excluding section 4) under the supervision of Notified body, **Szutest Teknik Kontrol ve Belgelendirme Hizmetleri Ticaret Limited Şirketi (2195).***

Notified Body : **Szutest Teknik Kontrol ve Belgelendirme Hizmetleri Ticaret Limited Şirketi**
Yukarı Dudullu Mahallesi Nato Yolu Caddesi Çam Sokak No: 7
Ümraniye/İstanbul, Turkey

EC Certificate No. : 2195-MED-1306401 (Issue date: 2013.04.25 / Valid date: 2018.03.04)

Place, Date of Issue : Republic of Korea / Feb. 12, 2013

Signature :



Yoo, Hag-Dong , **President**
On behalf of **HUBIT Co., Ltd.**

Attachment, Model/Type

a. PERFECT CLEAR

PERFECT CLEAR ROTH								
Maxillary	Torque	Ang .	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
Centrals	+11°	+5°	3.20	300-11	300-21	310-11	310-21	KOREA
Laterals	+8°	+9°	2.95	300-12	300-22	310-12	310-22	KOREA
Cuspids	-2°	+8°	3.50	300-13	300-23	310-13	310-23	KOREA
Hooks on Cuspids	-2°	+8°	3.50	300-13H	300-23H	310-13H	310-23H	KOREA
1 st Bicuspid	-7°	0°	3.20	300-14	300-24	310-14	310-24	KOREA
Hooks on 1 st Bicuspid	-7°	0°	3.20	300-14H	300-24H	310-14H	310-24H	KOREA
2 nd Bicuspid	-7°	0°	3.20	300-15	300-25	310-15	310-25	KOREA
Hooks on 2 nd Bicuspid	-7°	0°	3.20	300-15H	300-25H	310-15H	310-25H	KOREA
Mandibular								
Anteriors	0°	0°	2.85	300-41/42	300-31/32	310-41/42	310-31/32	KOREA
Cuspids	-11°	+7°	3.50	300-43	300-33	310-43	310-33	KOREA
Hooks on Cuspids	-11°	+7°	3.50	300-43H	300-33H	310-43H	310-33H	KOREA
1 st Bicuspid	-17°	+3°	3.20	300-44	300-34	310-44	310-34	KOREA
Hooks on 1 st Bicuspid	-17°	+3°	3.20	300-44H	300-34H	310-44H	310-34H	KOREA
1 st Bicuspid	-17°	0°	3.20	300-44A	300-34A	310-44A	310-34A	KOREA
Hooks on 1 st Bicuspid	-17°	0°	3.20	300-44HA	300-34HA	310-44HA	310-34HA	KOREA
2 nd Bicuspid	-21°	+6°	3.20	300-45	300-35	310-45	310-35	KOREA
Hooks on 2 nd Bicuspid	-21°	+6°	3.20	300-45H	300-35H	310-45H	310-35H	KOREA
2 nd Bicuspid	-21°	0°	3.20	300-45HA	300-35HA	310-45HA	310-35HA	KOREA
Hooks on 2 nd Bicuspid	-21°	0°	3.20	300-45A	300-35A	310-45A	310-35A	KOREA

ASSORTED KIT			.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)			300-030	310-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)			300-031	310-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)			300-300	310-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)			300-301	310-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)			300-330	310-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)			300-331	310-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)			300-050	310-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)			300-051	310-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)			300-053	310-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)			300-500	310-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)			300-501	310-501
Upper 5-5, Hooks on #3,4,5(11,12,13H,14H,15H 21,22,23H,24H,25H)			300-503	310-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)			300-550	310-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)			300-551	310-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H,31,32,33H,34H,35H,41,42,43H,44H,45H)			300-553	310-553

PERFECT CLEAR MBT								
Maxillary	Torque	Ang .	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
<i>Centrals</i>	+17°	+4°	3.20	320-11	320-21	330-11	330-21	KOREA
<i>Laterals</i>	+10°	+8°	2.95	320-12	320-22	330-12	330-22	KOREA
<i>Cuspids</i>	0°	+8°	3.50	320-13	320-23	330-13	330-23	KOREA
Hooks on Cuspids	0°	+8°	3.50	320-13H	320-23H	330-13H	330-23H	KOREA
<i>1st Bicuspid</i>	-7°	0°	3.20	320-14	320-24	330-14	330-24	KOREA
Hooks on 1st Bicuspid	-7°	0°	3.20	320-14H	320-24H	330-14H	330-24H	KOREA
<i>2nd Bicuspid</i>	-7°	0°	3.20	320-15	320-25	330-15	330-25	KOREA
Hooks on 2nd Bicuspid	-7°	0°	3.20	320-15H	320-25H	330-15H	330-25H	KOREA
Mandibular								
<i>Anteriors</i>	-6°	0°	2.85	320-41/42	320-31/32	330-41/42	330-31/32	KOREA
<i>Cuspids</i>	0°	+3°	3.50	320-43	320-33	330-43	330-33	KOREA
Hooks on Cuspids	0°	+3°	3.50	320-43H	320-33H	330-43H	330-33H	KOREA
<i>1st Bicuspid</i>	-12°	+2°	3.20	320-44	320-34	330-44	330-34	KOREA
Hooks on 1st Bicuspid	-12°	+2°	3.20	320-44H	320-34H	330-44H	330-34H	KOREA
<i>2nd Bicuspid</i>	-17°	+2°	3.20	320-45	320-35	330-45	330-35	KOREA
Hooks on 2nd Bicuspid	-17°	+2°	3.20	320-45H	320-35H	330-45H	330-35H	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	320-030	330-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	320-031	330-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	320-300	330-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	320-301	330-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	320-330	330-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	320-331	330-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	320-050	330-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	320-051	330-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)	320-053	330-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	320-500	330-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	320-501	330-501
Upper 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H 21,22,23H,24H,25H)	320-503	330-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	320-550	330-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	320-551	330-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H, 31,32,33H,34H,35H, 41,42,43H,44H,45H)	320-553	330-553

PERFECT CLEAR STANDARD								
Maxillary	Torque	Ang.	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
<i>Centrals</i>	0°	0°	3.20	340-11	340-21	350-11	350-21	KOREA
<i>Laterals</i>	0°	0°	3.00	340-12	340-22	350-12	350-22	KOREA
<i>Cuspids</i>	0°	0°	3.50	340-13	340-23	350-13	350-23	KOREA
Hooks on Cuspids	0°	0°	3.50	340-13H	340-23H	350-13H	350-23H	KOREA
<i>1st Bicuspid</i>	0°	0°	3.20	340-14	340-24	350-14	350-24	KOREA
<i>2nd Bicuspid</i>	0°	0°	3.20	340-15	340-25	350-15	350-25	KOREA
Mandibular								
<i>Anteriors</i>	0°	0°	2.85	340-41/42	340-31/32	350-41/42	350-31/32	KOREA
<i>Cuspids</i>	0°	0°	3.50	340-43	340-33	350-43	350-33	KOREA
Hooks on Cuspids	0°	0°	3.50	340-43H	340-33H	350-43H	350-33H	KOREA
<i>1st Bicuspid</i>	0°	0°	3.20	340-44	340-34	350-44	350-34	KOREA
<i>2nd Bicuspid</i>	0°	0°	3.20	340-45	340-35	350-45	350-35	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	340-030	350-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	340-031	350-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	340-300	350-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	340-301	350-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	340-330	350-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	340-331	350-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	340-050	350-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	340-051	350-051
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	340-500	350-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	340-501	350-501
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	340-550	350-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	340-551	350-551

b. PERFECT CLEAR II

PERFECT CLEAR II ROTH								
Maxillary	Torque	Angle	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
<i>Centrals</i>	+11°	+5°	3.20	200-11	200-21	210-11	210-21	KOREA
<i>Laterals</i>	+8°	+9°	2.95	200-12	200-22	210-12	210-22	KOREA
<i>Cuspids</i>	-2°	+8°	3.50	200-13	200-23	210-13	210-23	KOREA
<i>Hooks on Cuspids</i>	-2°	+8°	3.50	200-13H	200-23H	210-13H	210-23H	KOREA
<i>1st Bicuspid</i>	-7°	0°	3.20	200-14	200-24	210-14	210-24	KOREA
<i>Hooks on 1st Bicuspid</i>	-7°	0°	3.20	200-14H	200-24H	210-14H	210-24H	KOREA
<i>2nd Bicuspid</i>	-7°	0°	3.20	200-15	200-25	210-15	210-25	KOREA
<i>Hooks on 2nd Bicuspid</i>	-7°	0°	3.20	200-15H	200-25H	210-15H	210-25H	KOREA
Mandibular								
<i>Anteriors</i>	0°	0°	2.85	200-41/42	200-31/32	210-41/42	210-31/32	KOREA
<i>Cuspids</i>	-11°	+7°	3.50	200-43	200-33	210-43	210-33	KOREA
<i>Hooks on Cuspids</i>	-11°	+7°	3.50	200-43H	200-33H	210-43H	210-33H	KOREA
<i>1st Bicuspid</i>	-17°	+3°	3.20	200-44	200-34	210-44	210-34	KOREA
<i>Hooks on 1st Bicuspid</i>	-17°	+3°	3.20	200-44H	200-34H	210-44H	210-34H	KOREA
<i>1st Bicuspid</i>	-17°	0°	3.20	200-44A	200-34A	210-44A	210-34A	KOREA
<i>Hooks on 1st Bicuspid</i>	-17°	0°	3.20	200-44HA	200-34HA	210-44HA	210-34HA	KOREA
<i>2nd Bicuspid</i>	-21°	+6°	3.20	200-45	200-35	210-45	210-35	KOREA
<i>Hooks on 2nd Bicuspid</i>	-21°	+6°	3.20	200-45H	200-35H	210-45H	210-35H	KOREA
<i>2nd Bicuspid</i>	-21°	0°	3.20	200-45HA	200-35HA	210-45HA	210-35HA	KOREA
<i>Hooks on 2nd Bicuspid</i>	-21°	0°	3.20	200-45A	200-35A	210-45A	210-35A	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	200-030	210-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	200-031	210-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	200-300	210-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	200-301	210-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	200-330	210-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	200-331	210-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	200-050	210-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	200-051	210-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)	200-053	210-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	200-500	210-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	200-501	210-501
Upper 5-5, Hooks on #3,4,5(11,12,13H,14H,15H 21,22,23H,24H,25H)	200-503	210-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	200-550	210-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	200-551	210-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H,31,32,33H,34H,35H,41,42,43H,44H,45H)	200-553	210-553

PERFECT CLEAR II MBT								
Maxillary	Torque	Ang.	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
Centrals	+17°	+4°	3.20	220-11	220-21	230-11	230-21	KOREA
Laterals	+10°	+8°	2.95	220-12	220-22	230-12	230-22	KOREA
Cuspids	0°	+8°	3.50	220-13	220-23	230-13	230-23	KOREA
Hooks on Cuspids	0°	+8°	3.50	220-13H	220-23H	230-13H	230-23H	KOREA
1 st Bicuspid	-7°	0°	3.20	220-14	220-24	230-14	230-24	KOREA
Hooks on 1 st Bicuspid	-7°	0°	3.20	220-14H	220-24H	230-14H	230-24H	KOREA
2 nd Bicuspid	-7°	0°	3.20	220-15	220-25	230-15	230-25	KOREA
Hooks on 2 nd Bicuspid	-7°	0°	3.20	220-15H	220-25H	230-15H	230-25H	KOREA
Mandibular								
Anteriors	-6°	0°	2.85	220-41/42	220-31/32	230-41/42	230-31/32	KOREA
Cuspids	0°	+3°	3.50	220-43	220-33	230-43	230-33	KOREA
Hooks on Cuspids	0°	+3°	3.50	220-43H	220-33H	230-43H	230-33H	KOREA
1 st Bicuspid	-12°	+2°	3.20	220-44	220-34	230-44	230-34	KOREA
Hooks on 1 st Bicuspid	-12°	+2°	3.20	220-44H	220-34H	230-44H	230-34H	KOREA
2 nd Bicuspid	-17°	+2°	3.20	220-45	220-35	230-45	230-35	KOREA
Hooks on 2 nd Bicuspid	-17°	+2°	3.20	220-45H	220-35H	230-45H	230-35H	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	220-030	230-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	220-031	230-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	220-300	230-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	220-301	230-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	220-330	230-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	220-331	230-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	220-050	230-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	220-051	230-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)	220-053	230-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	220-500	230-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	220-501	230-501
Upper 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H 21,22,23H,24H,25H)	220-503	230-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	220-550	230-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	220-551	230-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H,31,32,33H,34H,35H,41,42,43H,44H,45H)	220-553	230-553

PERFECT CLEAR II STANDARD								
Maxillary	Torque	Ang.	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
Centrals	0°	0°	3.20	240-11	240-21	250-11	250-21	KOREA
Laterals	0°	0°	3.00	240-12	240-22	250-12	250-22	KOREA
Cuspids	0°	0°	3.50	240-13	240-23	250-13	250-23	KOREA
Hooks on Cuspids	0°	0°	3.50	240-13H	240-23H	250-13H	250-23H	KOREA
1 st Bicuspid	0°	0°	3.20	240-14	240-24	250-14	250-24	KOREA
2 nd Bicuspid	0°	0°	3.20	240-15	240-25	250-15	250-25	KOREA
Mandibular								
Anteriors	0°	0°	2.85	240-41/42	240-31/32	250-41/42	250-31/32	KOREA
Cuspids	0°	0°	3.50	240-43	240-33	250-43	250-33	KOREA
Hooks on Cuspids	0°	0°	3.50	240-43H	240-33H	250-43H	250-33H	KOREA
1 st Bicuspid	0°	0°	3.20	240-44	240-34	250-44	250-34	KOREA
2 nd Bicuspid	0°	0°	3.20	240-45	240-35	250-45	250-35	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	240-030	250-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	240-031	250-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	240-300	250-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	240-301	250-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	240-330	250-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	240-331	250-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	240-050	250-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	240-051	250-051
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	240-500	250-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	240-501	250-501
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	240-550	250-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	240-551	250-551

c. PERFECT CLEAR II (TRI-EDGE)

PERFECT CLEAR II (Tri-edge) ROTH						
Maxillary	Torque	Ang.	Width (mm)	.022 Slot		Approved for Sale in
				R	L	
Centrals	+11t	+51	3.20	210-11	210-21	KOREA
Laterals	+8t	+9t	2.95	210-12	210-22	KOREA
Hooks on Cuspids	-2s	+8s	3.50	210-13H	210-23H	KOREA
1 st Bicuspids	-7i	07	3.20	210-14	210-24	KOREA
Hooks on 1 st Bicuspids	-7i	07	3.20	210-14H	210-24H	KOREA
2 nd Bicuspids	-7i	07	3.20	210-15	210-25	KOREA
Hooks on 2 nd Bicuspids	-7i	07	3.20	210-15H	210-25H	KOREA
Mandibular						
Anteriors	0t	0t	2.85	210-41/42	210-31/32	KOREA
Hooks on Cuspids	-11p	+71	3.50	210-43HT	210-33HT	KOREA
1 st Bicuspids	-17c	+37	3.20	210-44T	210-34T	KOREA
2 nd Bicuspids	-21c	+61	3.20	210-45T	210-35T	KOREA

ASSORTED KIT	ROTH 022	MBT 022
Upper/Lower 3-3, Hooks on Cuspids	210-331T	230-331T
Upper/Lower 5-5, Hooks on Cuspids	210-551T	230-551T

PERFECT CLEAR II (Tri-edge) MBT						
Maxillary	Torque	Ang.	Width (mm)	.022 Slot		Approved for Sale in
				R	L	
Centrals	+17°	+4°	3.20	230-11	230-21	KOREA
Laterals	+10°	+8°	2.95	230-12	230-22	KOREA
Hooks on Cuspids	0°	+8°	3.50	230-13H	230-23H	KOREA
1 st Bicuspids	-7°	0°	3.20	230-14	230-24	KOREA
Hooks on 1 st Bicuspids	-7°	0°	3.20	230-14H	230-24H	KOREA
2 nd Bicuspids	-7°	0°	3.20	230-15	230-25	KOREA
Hooks on 2 nd Bicuspids	-7°	0°	3.20	230-15H	230-25H	KOREA
Mandibular						
Anteriors	-6°	0°	2.85	230-41/42	230-31/32	KOREA
Hooks on Cuspids	0°	+3°	3.50	230-43HT	230-33HT	KOREA
1 st Bicuspids	-12°	+2°	3.20	230-44T	230-34T	KOREA
2 nd Bicuspids	-17°	+2°	3.20	230-45T	230-35T	KOREA

d. WOW Bracket

wOW ROTH								
Maxillary	Torque	Ang.	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
<i>Centrals</i>	+11°	+5°	3.20	102-11	102-21	112-11	112-21	KOREA
<i>Laterals</i>	+8°	+9°	2.95	102-12	102-22	112-12	112-22	KOREA
<i>Cuspids</i>	-2°	+8°	3.50	102-13(H)	102-23(H)	112-13(H)	112-23(H)	KOREA
<i>1st Bicuspids</i>	-7°	0°	3.20	102-14(H)	102-24(H)	112-14(H)	112-24(H)	KOREA
<i>2nd Bicuspids</i>	-7°	0°	3.20	102-15(H)	102-25(H)	112-15(H)	112-25(H)	KOREA
Mandibular								
<i>Anteriors</i>	0°	0°	2.85	102-41/42	102-31/32	112-41/42	112-31/32	KOREA
<i>Cuspids</i>	-11°	+7°	3.50	102-43(H)	102-33(H)	112-43(H)	112-33(H)	KOREA
<i>1st Bicuspids</i>	-17°	0°	3.20	102-44(H)	102-34(H)	112-44(H)	112-34(H)	KOREA
<i>2nd Bicuspids</i>	-17°	0°	3.20	102-45(H)	102-35(H)	112-45(H)	112-35(H)	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	102-030	112-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	102-031	112-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	102-300	112-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	102-301	112-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	102-330	112-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	102-331	112-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	102-050	112-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	102-051	112-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)	102-053	112-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	102-500	112-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	102-501	112-501
Upper 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H 21,22,23H,24H,25H)	102-503	112-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	102-550	112-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	102-551	112-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H, 31,32,33H,34H,35H, 41,42,43H,44H,45H)	102-553	112-553

wow MBT								
<i>Maxillary</i>	Torque	Ang.	Width (mm)	.018 Slot		.022 Slot		Approved for Sale in
				R	L	R	L	
<i>Centrals</i>	+17°	+4°	3.20	122-11	122-21	132-11	132-21	KOREA
<i>Laterals</i>	+10°	+8°	2.95	122-12	122-22	132-12	132-22	KOREA
<i>Cuspids</i>	0°	+8°	3.50	122-13(H)	122-23(H)	132-13(H)	132-23(H)	KOREA
<i>1st Bicuspid</i>	-7°	0°	3.20	122-14(H)	122-24(H)	132-14(H)	132-24(H)	KOREA
<i>2nd Bicuspid</i>	-7°	0°	3.20	122-15(H)	122-25(H)	132-15(H)	132-25(H)	KOREA
<i>Mandibular</i>								
<i>Anteriors</i>	-6°	0°	2.85	122-41/42	122-31/32	132-41/42	132-31/32	KOREA
<i>Cuspids</i>	0°	+3°	3.50	122-43(H)	122-33(H)	132-43(H)	132-33(H)	KOREA
<i>1st Bicuspid</i>	-12°	0°	3.20	122-44(H)	122-34(H)	132-44(H)	132-34(H)	KOREA
<i>2nd Bicuspid</i>	-12°	0°	3.20	122-45(H)	122-35(H)	132-45(H)	132-35(H)	KOREA

ASSORTED KIT	.018 Slot	.022 Slot
Lower 3-3, Without Hooks on Cuspids(31,32,33, 41,42,43)	122-030	132-030
Lower 3-3, Hooks on Cuspids(31,32,33H, 41,42,43H)	122-031	132-031
Upper 3-3, Without Hooks on Cuspids(11,12,13, 21,22,23)	122-300	132-300
Upper 3-3, Hooks on Cuspids(11,12,13H, 21,22,23H)	122-301	132-301
Upper and Lower 3-3, Without Hooks on Cuspids (11,12,13, 21,22,23, 31,32,33, 41,42,43)	122-330	132-330
Upper and Lower 3-3, Hooks on Cuspids (11,12,13H, 21,22,23H, 31,32,33H, 41,42,43H)	122-331	132-331
Lower 5-5, Without Hooks on Cuspids(31,32,33,34,35 41,42,43,44,45)	122-050	132-050
Lower 5-5, Hooks on Cuspids(31,32,33H,34,35 41,42,43H,44,45)	122-051	132-051
Lower 5-5, Hooks on #3,4,5 (31,32,33H,34H,35H 41,42,43H,44H,45H)	122-053	132-053
Upper 5-5, Without Hooks on Cuspids(11,12,13,14,15 21,22,23,24,25)	122-500	132-500
Upper 5-5, Hooks on Cuspids(11,12,13H,14,15 21,22,23H,24,25)	122-501	132-501
Upper 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H 21,22,23H,24H,25H)	122-503	132-503
Upper and Lower 5-5, Without Hooks on Cuspids (11,12,13,14,15 21,22,23,24,25, 31,32,33,34,35 41,42,43,44,45)	122-550	132-550
Upper and Lower 5-5, Hooks on Cuspids (11,12,13H,14,15 21,22,23H,24,25, 31,32,33H,34,35 41,42,43H,44,45)	122-551	132-551
Upper and Lower 5-5, Hooks on #3,4,5 (11,12,13H,14H,15H,21,22,23H,24H,25H,31,32,33H,34H,35H, 41,42,43H,44H,45H)	122-553	132-553

End of File.